


Managing a Workspace

Open Workspace




The **Open Workspace...** menu item is only enabled if an LMD file is open in the Workspace. If no LMD file is open, the **Open Workspace...** button is disabled.

To open an existing Workspace file click the **Open Workspace...** option (or press the **Ctrl+F** keys).

The Windows **Open** dialog is displayed, listing the files present in the last used folder or **My Documents** if this is the first time you have used the software.

Select the desired workspace file and click the **Open** button to load the workspace file with the currently open LMD file.

If the current folder is not the one containing your workspace files, you need to browse to the correct folder.



The files of type displays **Workspace Files (*.xmlw)** and can be changed to display **All Files (*.*)**.


Double clicking on a file in this dialog has the same functionality as selecting it and clicking the **Open** button.

If the parameters in the selected Workspace do not match those of the current LMD the Parameter Resolver dialog is displayed allowing you to select parameters matching the current workspace.

Save Workspace

A VenturiOne workspace includes plots, regions, gates, color precedence, selected statistics and displayed parameters. A workspace can be saved for later opening.

To save a current workspace with the current workspace file name, click the **Save Workspace** option or press the Ctrl+S keys:




The file is always saved with an **xmlw** extension.


If a modified untitled workspace file is to be saved, a **Save As** dialog is displayed, allowing you to save the current workspace with a new name.

Save Workspace As...

To save the current workspace with a new file name, click the **Save Workspace As...** option (or press the Ctrl + D keys)



The **Save As** dialog is displayed to allow the assignment of a relevant file name as illustrated below:



The dialog defaults to the last used folder or **My Documents** if there is no last used folder information.

The **File name:** field contains **workspace** by default.

Workspace Files (*.xmlw) appears in the **Save as type:** drop-down menu.

The **Save as type:** drop-down menu can also be changed to '**All Files**' (*.*) .

You can enter a suitable name for the workspace.

Click the **Save** button to save the current workspace with the chosen name. If the selected file name already exists you will have the option to overwrite or change the file name.

Click the **Cancel** button or **Close** button (X top right) on the **Save As** dialog, to close the dialog without saving the workspace file.